





BLK Antibody

Product Code	CSB-PA002714ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51451
Immunogen	Recombinant Human Tyrosine-protein kinase Blk protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Blk (EC 2.7.10.2) (B lymphocyte kinase) (p55-Blk), BLK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BLK

Image



Western blot
 All lanes: Tyrosine-protein kinase Blk antibody at 10µg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58 kDa
 Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA002714ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.