




CUL4B Antibody

Product Code	CSB-PA615708ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13620
Immunogen	Recombinant Human Cullin-4B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cullin-4B (CUL-4B), CUL4B, KIAA0695
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CUL4B

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA615708ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.