




ACAA1 Antibody

Product Code	CSB-PA001117ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09110
Immunogen	Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal protein (27-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-ketoacyl-CoA thiolase, peroxisomal (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase) (Peroxisomal 3-oxoacyl-CoA thiolase), ACAA1, ACAA PTHIO
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACAA1

Image



Western blot

All lanes: 3-ketoacyl-CoA thiolase, peroxisomal antibody at 8µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 45, 35 kDa
 Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.