




CALHM1 Antibody

Product Code	CSB-PA808532ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU99
Immunogen	Recombinant Human Calcium homeostasis modulator protein 1 protein (202-346AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium homeostasis modulator protein 1 (Protein FAM26C), CALHM1, FAM26C
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CALHM1

Image



Western blot
 All lanes: Calcium homeostasis modulator protein 1 antibody at 7µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 38 kDa
 Observed band size: 38 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.