





EME1 Antibody

Product Code	CSB-PA856897LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AY2
Immunogen	Recombinant Human Crossover junction endonuclease EME1 protein (1-583AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Crossover junction endonuclease EME1 (EC 3.1.22.-) (MMS4 homolog) (hMMS4), EME1, MMS4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EME1

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA856897LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA856897LA01HU at dilution of 1:100