



Pectate lyase 5 Antibody, Biotin conjugated

Product Code	CSB-PA329754HD01BYC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27759
Immunogen	Recombinant Ambrosia artemisiifolia Pectate lyase 5 protein (26-396AA)
Raised In	Rabbit
Species Reactivity	Ambrosia artemisiifolia
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 5 (EC 4.2.2.2) (Antigen Amb a I) (Antigen E) (AgE) (Pollen allergen Amb a 1.1) (allergen Amb a 1.1)
Immunogen Species	Ambrosia artemisiifolia (Short ragweed)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.