



Botulinum neurotoxin type C1 Antibody, Biotin conjugated

Product Code	CSB-PA325596LD01CLQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18640
Immunogen	Recombinant Clostridium botulinum Botulinum neurotoxin type C1 protein (2-428AA)
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Botulinum neurotoxin type C1 (BoNT/C1) (EC 3.4.24.69) (Bontoxilysin-C1) [Cleaved into: Botulinum neurotoxin C1 light chain; Botulinum neurotoxin C1 heavy chain]
Immunogen Species	Clostridium botulinum
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.