



# stxB2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA357663HD01ECC
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09386
<b>Immunogen</b>	Recombinant Enterobacteria phage 933W Shiga-like toxin 2 subunit B protein (20-89AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Enterobacteria phage 933W
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Shiga-like toxin 2 subunit B (SLT-2 B subunit) (SLT-2b) (SLT-IIb) (Verocytotoxin 2 subunit B) (Verotoxin 2 subunit B), stxB2, stx2B
<b>Immunogen Species</b>	Enterobacteria phage 933W (Bacteriophage 933W)
<b>Research Area</b>	Others
<b>Target Names</b>	stxB2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.