



stxB Antibody, HRP conjugated

Product Code	CSB-PA300809HB01EKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69179
Immunogen	Recombinant Enterobacteria phage H19B Shiga-like toxin 1 subunit B protein (21-89AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage H19B
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Shiga-like toxin 1 subunit B (SLT-1 B subunit) (SLT-1b) (SLT-1b) (Verocytotoxin 1 subunit B) (Verotoxin 1 subunit B), stxB, sltB
Immunogen Species	Enterobacteria phage H19B (Bacteriophage H19B)
Research Area	Others
Target Names	stxB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.