





TFEB Antibody

Product Code	CSB-PA023429ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19484
Immunogen	Recombinant Human Transcription factor EB protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TFEB

Image



Western blot
 All lanes: Transcription factor EB antibody at 3µg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: Mouse brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 53, 44 kDa
 Observed band size: 53 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA023429ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.