




SLC27A4 Antibody

Product Code	CSB-PA761274ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1M0
Immunogen	Recombinant Human Long-chain fatty acid transport protein 4 protein (494-643AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain fatty acid transport protein 4 (FATP-4) (Fatty acid transport protein 4) (EC 6.2.1.-) (Solute carrier family 27 member 4), SLC27A4, ACSVL4 FATP4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC27A4

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA761274ESR2HU at dilution of 1:100