





GRB7 Antibody

Product Code	CSB-PA009890ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14451
Immunogen	Recombinant Human Growth factor receptor-bound protein 7 protein (233-532AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth factor receptor-bound protein 7 (B47) (Epidermal growth factor receptor GRB-7) (GRB7 adapter protein), GRB7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRB7

Image



Western blot
 All lanes: Growth factor receptor-bound protein 7 antibody at 2µg/ml + Rat brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 60, 50, 63, 56 kDa
 Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human esophagus tissue using CSB-PA009890ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.