





GRB7 Antibody

Product Code	CSB-PA009890ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14451
Immunogen	Recombinant Human Growth factor receptor-bound protein 7 protein (233-532AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth factor receptor-bound protein 7 (B47) (Epidermal growth factor receptor GRB-7) (GRB7 adapter protein), GRB7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRB7

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA009890ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human esophagus tissue using CSB-PA009890ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.