




PHYH Antibody

Product Code	CSB-PA017946DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14832
Immunogen	Recombinant Human Phytanoyl-CoA dioxygenase, peroxisomal protein (1-338AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phytanoyl-CoA dioxygenase, peroxisomal (EC 1.14.11.18) (Phytanic acid oxidase) (Phytanoyl-CoA alpha-hydroxylase) (PhyH), PHYH, PAHX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PHYH

Image



Western blot


All lanes: Phytanoyl-CoA dioxygenase, peroxisomal antibody at 4µg/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 39, 28 kDa

Observed band size: 39 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA017946DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.