





PRND Antibody

Product Code	CSB-PA890741ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKY0
Immunogen	Recombinant Human Prion-like protein doppel protein (20-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prion-like protein doppel (PrPLP) (Prion protein 2), PRND, DPL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRND

Image



Western blot
 All lanes: Prion-like protein doppel antibody at 2µg/ml + Mouse spleen tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 20 kDa
 Observed band size: 20 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA890741ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.