




PRMT1 Antibody


Product Code	CSB-PA859113ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99873
Immunogen	Recombinant Human Protein arginine N-methyltransferase 1 protein (20-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein arginine N-methyltransferase 1 (EC 2.1.1.319) (Histone-arginine N-methyltransferase PRMT1) (Interferon receptor 1-bound protein 4), PRMT1, HMT2 HRMT1L2 IR1B4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRMT1

Image



Western blot

All lanes: Protein arginine N-methyltransferase 1 antibody at 2µg/ml + NIH/3T3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 42, 40, 41 kDa
 Observed band size: 42 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA859113ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.