





SMYD3 Antibody

Product Code	CSB-PA884482ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7B4
Immunogen	Recombinant Human Histone-lysine N-methyltransferase SMYD3 protein (199-428AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone-lysine N-methyltransferase SMYD3 (EC 2.1.1.43) (SET and MYND domain-containing protein 3) (Zinc finger MYND domain-containing protein 1), SMYD3, ZMYND1 ZNFN3A1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SMYD3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA884482ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA884482ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.