



Trim21 Antibody, HRP conjugated

Product Code	CSB-PA733816LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62191
Immunogen	Recombinant Mouse E3 ubiquitin-protein ligase TRIM21 protein (1-470AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM21 (EC 2.3.2.27) (52 kDa Ro protein) (52 kDa ribonucleoprotein autoantigen Ro/SS-A) (RING-type E3 ubiquitin transferase TRIM21) (Ro(SS-A)) (Sjogren syndrome type A antigen) (SS-A) (Tripartite motif-containing protein 21), Trim21, Ro52 Ssa1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Trim21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.