



# Trim21 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA733816LD01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q62191
<b>Immunogen</b>	Recombinant Mouse E3 ubiquitin-protein ligase TRIM21 protein (1-470AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	E3 ubiquitin-protein ligase TRIM21 (EC 2.3.2.27) (52 kDa Ro protein) (52 kDa ribonucleoprotein autoantigen Ro/SS-A) (RING-type E3 ubiquitin transferase TRIM21) (Ro(SS-A)) (Sjogren syndrome type A antigen) (SS-A) (Tripartite motif-containing protein 21), Trim21, Ro52 Ssa1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Others
<b>Target Names</b>	Trim21
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.