




SOAT1 Antibody

Product Code	CSB-PA022385DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35610
Immunogen	Recombinant Human Sterol O-acyltransferase 1 protein (1-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sterol O-acyltransferase 1 (EC 2.3.1.26) (Acyl-coenzyme A:cholesterol acyltransferase 1) (ACAT-1) (Cholesterol acyltransferase 1), SOAT1, ACACT ACACT1 ACAT ACAT1 SOAT STAT
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SOAT1

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA022385DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.