



H671_3g9701 Antibody, HRP conjugated

Product Code	CSB-PA006187LB01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G3I4W7
Immunogen	Recombinant Cricetulus griseus Cathepsin D protein (65-408AA)
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin D, H671_3g9701
Immunogen Species	Cricetulus griseus
Research Area	Others
Target Names	H671_3g9701
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.