



tetX Antibody, Biotin conjugated

Product Code	CSB-PA356416LD01CMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04958
Immunogen	Recombinant Clostridium tetani Tetanus toxin protein (123-573AA)
Raised In	Rabbit
Species Reactivity	Clostridium tetani
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetanus toxin (EC 3.4.24.68) (Tentoxylisin) [Cleaved into: Tetanus toxin light chain (Tetanus toxin chain L); Tetanus toxin heavy chain (Tetanus toxin chain H)], tetX
Immunogen Species	Clostridium tetani (strain Massachusetts / E88)
Research Area	Others
Target Names	tetX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.