



cssB Antibody, HRP conjugated

Product Code	CSB-PA347417HB01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53510
Immunogen	Recombinant Escherichia coli CS6 fimbrial subunit B protein (22-167AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CS6 fimbrial subunit B, cssB
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	cssB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.