



# tox B Antibody, HRP conjugated

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA321589HB01DUN                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | P18177                                                                           |
| Immunogen           | Recombinant Peptoclostridium difficile Toxin B protein (1751-2366AA)             |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Peptoclostridium difficile                                                       |
| Tested Applications | ELISA                                                                            |
| Form                | Liquid                                                                           |
| Conjugate           | HRP                                                                              |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified                                                         |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Alias               | Toxin B (EC 3.4.22.-), tox B, tcdB                                               |
| Immunogen Species   | Peptoclostridium difficile                                                       |
| Research Area       | Others                                                                           |
| Target Names        | tox B                                                                            |