




DCAF7 Antibody

Product Code	CSB-PA026056ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61962
Immunogen	Recombinant Human DDB1- and CUL4-associated factor 7 protein (1-342AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DDB1- and CUL4-associated factor 7 (WD repeat-containing protein 68) (WD repeat-containing protein An11 homolog), DCAF7, HAN11 WDR68
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DCAF7

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA026056ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.