




SNAP25 Antibody

Product Code	CSB-PA021873ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60880
Immunogen	Recombinant Human Synaptosomal-associated protein 25 protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptosomal-associated protein 25 (SNAP-25) (Super protein) (SUP) (Synaptosomal-associated 25 kDa protein), SNAP25, SNAP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNAP25

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA021873ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.