




GCK Antibody

Product Code	CSB-PA009319ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35557
Immunogen	Recombinant Human Glucokinase protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucokinase (EC 2.7.1.2) (Hexokinase type IV) (HK IV) (Hexokinase-4) (HK4) (Hexokinase-D), GCK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GCK

Image



Western blot


All lanes: Glucokinase antibody at 8µg/ml + Jurkat whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA009319ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA009319ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.