




ALPP Antibody

Product Code	CSB-PA001632ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05187
Immunogen	Recombinant Human Alkaline phosphatase, placental type protein (260-520AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkaline phosphatase, placental type (EC 3.1.3.1) (Alkaline phosphatase Regan isozyme) (Placental alkaline phosphatase 1) (PLAP-1), ALPP, PLAP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ALPP

Image



Western blot

All lanes: Alkaline phosphatase, placental type antibody at 8µg/ml

Lane 1: A549 whole cell lysate


Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA001632ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.