




SLC25A13 Antibody

Product Code	CSB-PA892444ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJS0
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein Aralar2 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein Aralar2 (Citrin) (Mitochondrial aspartate glutamate carrier 2) (Solute carrier family 25 member 13), SLC25A13, ARALAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC25A13

Image



Western blot

All lanes: Calcium-binding mitochondrial carrier protein Aralar2 antibody at 4µg/ml

Lane 1: Mouse liver tissue

Lane 2: Mouse kidney tissue


Lane 2: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 74 kDa

Observed band size: 74 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA892444ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.