



Tnnt2 Antibody

Product Code	CSB-PA024016HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50752
Immunogen	Recombinant Mouse Troponin T, cardiac muscle protein (2-301AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin T, cardiac muscle (TnTc) (Cardiac muscle troponin T) (cTnT), Tnnt2
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Tnnt2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.