




NCS1 Antibody

Product Code	CSB-PA008983ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62166
Immunogen	Recombinant Human Neuronal calcium sensor 1 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal calcium sensor 1 (NCS-1) (Frequenin homolog) (Frequenin-like protein) (Frequenin-like ubiquitous protein), NCS1, FLUP FREQ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NCS1

Image



Western blot

All lanes: Neuronal calcium sensor 1 antibody at 10μg/ml

Lane 1: Hela whole cell lysate

Lane 2: Rat kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 22, 20 kDa

Observed band size: 22 kDa