




# NCS1 Antibody

<b>Product Code</b>	CSB-PA008983ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P62166
<b>Immunogen</b>	Recombinant Human Neuronal calcium sensor 1 protein (1-190AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Neuronal calcium sensor 1 (NCS-1) (Frequenin homolog) (Frequenin-like protein) (Frequenin-like ubiquitous protein), NCS1, FLUP FREQ
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	NCS1

## Image



### Western blot

All lanes: Neuronal calcium sensor 1 antibody at 10µg/ml

Lane 1: HeLa whole cell lysate

Lane 2: Rat kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 22, 20 kDa

Observed band size: 22 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.