



# gag Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA15019C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04591
<b>Immunogen</b>	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag polyprotein protein (2-132AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human immunodeficiency virus type 1 group M subtype B
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Gag polyprotein (Pr55Gag) [Cleaved into: Matrix protein p17 (MA); Capsid protein p24 (CA); Spacer peptide 1 (SP1) (p2); Nucleocapsid protein p7 (NC); Spacer peptide 2 (SP2) (p1); p6-gag], gag
<b>Immunogen Species</b>	Human immunodeficiency virus type 1 group M subtype B
<b>Research Area</b>	Others
<b>Target Names</b>	gag
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.