



gag Antibody, Biotin conjugated

Product Code	CSB-PA15019D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04591
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag polyprotein protein (2-132AA)
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus type 1 group M subtype B
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gag polyprotein (Pr55Gag) [Cleaved into: Matrix protein p17 (MA); Capsid protein p24 (CA); Spacer peptide 1 (SP1) (p2); Nucleocapsid protein p7 (NC); Spacer peptide 2 (SP2) (p1); p6-gag], gag
Immunogen Species	Human immunodeficiency virus type 1 group M subtype B
Research Area	Others
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.