



TSHB Antibody

Product Code	CSB-PA11769A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01222
Immunogen	Recombinant Human Thyrotropin subunit beta protein (21-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyrotropin subunit beta (Thyroid-stimulating hormone subunit beta) (TSH-B) (TSH-beta) (Thyrotropin beta chain) (Thyrotropin alfa), TSHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSHB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.