



TSHB Antibody, FITC conjugated

Product Code	CSB-PA11769C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01222
Immunogen	Recombinant Human Thyrotropin subunit beta protein (21-132AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyrotropin subunit beta (Thyroid-stimulating hormone subunit beta) (TSH-B) (TSH-beta) (Thyrotropin beta chain) (Thyrotropin alfa), TSHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSHB