



PSMD10 Antibody, HRP conjugated

Product Code	CSB-PA018899LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75832
Immunogen	Recombinant Human 26S proteasome non-ATPase regulatory subunit 10 protein (1-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 10 (26S proteasome regulatory subunit p28) (Gankyrin) (p28(GANK)), PSMD10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.