




BTN2A2 Antibody

Product Code	CSB-PA848825LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVV5
Immunogen	Recombinant Human Butyrophilin subfamily 2 member A2 protein (33-262AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Butyrophilin subfamily 2 member A2, BTN2A2, BT2.2 BTF2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BTN2A2

Image



Western blot

All lanes: BTN2A2 antibody at 10µg/ml

Lane 1: MCF-7 whole cell lysate


Lane 2: NIH/3T3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60, 38, 29, 46, 36 kDa

Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded

human small intestine tissue using CSB-PA848825LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.