



FZR1 Antibody, Biotin conjugated

Product Code	CSB-PA892473LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM11
Immunogen	Recombinant Human Fizzy-related protein homolog protein (1-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fizzy-related protein homolog (Fzr) (CDC20-like protein 1) (Cdh1/Hct1 homolog) (hCDH1), FZR1, CDH1 FYR FZR KIAA1242
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FZR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.