



GCY1 Antibody, HRP conjugated

Product Code	CSB-PA321197HB01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14065
Immunogen	Recombinant Saccharomyces cerevisiae Glycerol 2-dehydrogenase (NADP(+)) protein (1-312AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol 2-dehydrogenase (NADP(+)) (EC 1.1.1.156) (Galactose-inducible crystallin-like protein 1), GCY1, GCY
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	GCY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.