



Polygalacturonase Antibody, FITC conjugated

Product Code	CSB-PA337345LC01DYO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43212
Immunogen	Recombinant Cryptomeria japonica Polygalacturonase protein (46-433AA)
Raised In	Rabbit
Species Reactivity	Cryptomeria japonica
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polygalacturonase (PG) (EC 3.2.1.15) (Allergen Cry j II) (Major pollen allergen Cry j 2) (Pectinase) (allergen Cry j 2)
Immunogen Species	Cryptomeria japonica (Japanese cedar) (Cupressus japonica)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.