



G Antibody, Biotin conjugated

Product Code	CSB-PA365937LD01HPO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03423
Immunogen	Recombinant Human respiratory syncytial virus A Major surface glycoprotein G protein (67-298AA)
Raised In	Rabbit
Species Reactivity	Human respiRatory syncytial virus A
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major surface glycoprotein G (Attachment glycoprotein G) (Membrane-bound glycoprotein) (mG), G
Immunogen Species	Human respiratory syncytial virus A (strain A2)
Research Area	Others
Target Names	G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.