



wnt8 Antibody

| | |
|---------------------|--|
| Product Code | CSB-PA340936LA01XBE |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P28026 |
| Immunogen | Recombinant Xenopus laevis Protein Wnt-8 protein (23-358AA) |
| Raised In | Rabbit |
| Species Reactivity | Xenopus laevis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein Wnt-8 (XWnt-8), wnt8 |
| Immunogen Species | Xenopus laevis (African clawed frog) |
| Research Area | Others |
| Target Names | wnt8 |