



wnt8 Antibody

Product Code	CSB-PA340936LA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28026
Immunogen	Recombinant Xenopus laevis Protein Wnt-8 protein (23-358AA)
Raised In	Rabbit
Species Reactivity	Xenopus laevis
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wnt-8 (XWnt-8), wnt8
Immunogen Species	Xenopus laevis (African clawed frog)
Research Area	Others
Target Names	wnt8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.