





RABGEF1 Antibody

Product Code	CSB-PA883453LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJ41
Immunogen	Recombinant Human Rab5 GDP/GTP exchange factor protein (1-491AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab5 GDP/GTP exchange factor (RAP1) (Rabaptin-5-associated exchange factor for Rab5) (Rabex-5), RABGEF1, RABEX5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RABGEF1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA883453LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA883453LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.