



# AC1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA324434LC01ADX
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14982
<b>Immunogen</b>	Recombinant African cassava mosaic virus Replication-associated protein (1-358AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	African cassava mosaic virus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Replication-associated protein (Rep) (EC 2.7.7.-) (EC 3.1.21.-) (40.4 kDa protein) (Protein AC1) (Protein AL1)
<b>Immunogen Species</b>	African cassava mosaic virus (isolate West Kenyan 844) (ACMV) (Cassava latent virus (isolate West Kenyan 844))
<b>Research Area</b>	Others
<b>Target Names</b>	AC1