



AC1 Antibody, Biotin conjugated

Product Code	CSB-PA324434LD01ADX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14982
Immunogen	Recombinant African cassava mosaic virus Replication-associated protein (1-358AA)
Raised In	Rabbit
Species Reactivity	African cassava mosaic virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Replication-associated protein (Rep) (EC 2.7.7.-) (EC 3.1.21.-) (40.4 kDa protein) (Protein AC1) (Protein AL1)
Immunogen Species	African cassava mosaic virus (isolate West Kenyan 844) (ACMV) (Cassava latent virus (isolate West Kenyan 844))
Research Area	Others
Target Names	AC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.