




DYNLRB1 Antibody

Product Code	CSB-PA865093HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP97
Immunogen	Recombinant Human Dynein light chain roadblock-type 1 protein (3-96AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein light chain roadblock-type 1 (Bithoraxoid-like protein) (BLP) (Dynein light chain 2A, cytoplasmic) (Dynein-associated protein Km23) (Roadblock domain-containing protein 1), DYNLRB1, BITH DNCL2A DNLC2A ROBLD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DYNLRB1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA865093HA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.