



DYNLRB1 Antibody, FITC conjugated

Product Code	CSB-PA865093HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP97
Immunogen	Recombinant Human Dynein light chain roadblock-type 1 protein (3-96AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein light chain roadblock-type 1 (Bithoraxoid-like protein) (BLP) (Dynein light chain 2A, cytoplasmic) (Dynein-associated protein Km23) (Roadblock domain-containing protein 1), DYNLRB1, BITH DNCL2A DNLC2A ROBLD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DYNLRB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.