



GSN Antibody, HRP conjugated

Product Code	CSB-PA009965HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06396
Immunogen	Recombinant Human Gelsolin protein (28-782AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gelsolin (AGEL) (Actin-depolymerizing factor) (ADF) (Brevin), GSN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.